

SonoTip® TopGain

DE	Biopsie-Nadel-System für die endobronchial sono-graphisch geführte transbronchiale Feinnadel-Biopsie (EBUS-TBNB) für Olympus und Fujifilm EBUS Bronchoskope	Gebrauchs-anweisung	6
GB	Biopsy needle system for Endobronchial Ultrasound-Guided Trans-bronchial Fine Needle Biopsy (EBUS-TBNB) for Olympus and Fujifilm EBUS Bronchoscopes	Instruction Manual	11
FR	Système Aiguille d'endoscopie ultrasonique pour biopsie par aiguille fine transbronchiale et endobronchiale (EBUS-TBNB) guidée par ultrason pour bronchoscope Olympus et Fujifilm EBUS	Manuel d'utilisation	16

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Warning !

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Intended Purpose

The SonoTip® TopGain EBUS Needle System is intended for Ultrasonically Guided Fine Needle Biopsy, (FNB) of submucosal and extraluminal lesions of the Tracheobronchial Tree and Gastrointestinal Tract (e.g., lymph nodes, abnormal tissue in the mediastinum).

! The presence of a cytologist during extraluminal EBUS guided TBNB may reduce the number of fine needle aspiration passes required from each lesion.

Contraindications

There are no special contraindications existing for this product.

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the targeted site.

Warning !

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury to the patient and/or user. It is essential to follow the instructions enclosed, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization.

! Possible complications from the extraluminal endobronchial ultrasound-guided TBNB may include but are not limited to infection, hemorrhage, perforation and tumor seeding.

! Failure to follow these instructions may result in instrument malfunction and/or damage!

Storage

Please store the instrument at room temperature, keep dry and avoid exposure to sunlight.

General Safety Instruction

- Wear protective garments (i.e., gloves, protective eyewear, gowns, etc.) to assist in the protection against cross-contamination and injury.
- Do not use this instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. The working channel of the endoscope must be at least "0.2 mm" larger than the SonoTip® TopGain EBUS outer diameter.
- ! Never use instruments which do not comply with the technical specifications.
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.

Inspecting and Testing

- Inspect the outer packaging for any defects or damaged seals.
- ! If damaged or an open seal is detected, sterilization may be compromised and the needle system should be discarded and replaced.
- Check the date of sterilization.
- ! Never use products that indicate an expired sterilization date.


Preparation of the Aspiration Syringe

- Open package and remove the sterile aspiration syringe.
- Check to ensure that the stopcock valve is securely attached to the syringe.
- Ensure syringe plunger is fully encased in the syringe chamber.
- Turn the stopcock valve to its closed position.
- Pull the plunger completely back. Turn the plunger clockwise until the plunger is locked, thereby creating negative suction pressure.
- Set the prepared syringe aside until you are ready to use it for biopsy.

Preparing the SonoTip® TopGain EBUS needle


- Remove the SonoTip® TopGain EBUS needle from its packaging and protective sleeve by sliding off the plastic sleeve.
- Inspect the entire system to insure that there are no defects and that the device functions properly. Pay particular attention to the following:
 - a) Verify that the needle **(10)** and the plastic sheath **(9)** are free of any noticeable kinks and/or defects.
 - b) Insure that the connection between plastic sheath **(9)** and handle body is secure.

Verify that the Sliding Sheath Handle **(7)** is fully retracted. **Picture 1** shows the correct position on the Scale for proper Sheath Adjustment **(5)**.

- ! If it is determined that the Sliding Sheath Handle **(7)** is not fixed in the correct position as shown in **Picture 1**, please unlock the blue Sheath Adjuster **(6)** by turning it counter-clockwise. Retract the Sliding Sheath Handle **(7)** until it stops as illustrated in **Picture 1**. Turn the Sheath Adjuster clockwise to lock the Sliding Sheath Handle in position.

Ensure that the Needle **(10)** and Nitinol Stylet are securely retracted into the Plastic Sheath **(9)** before advancing the instrument into the accessory channel of the endoscope. **Picture 2** shows the correct position on the Scale Needle Adjustment **(4)**.


- ! If it is determined that the Scale for Adjustment **(4)** is not in the correct position as shown in **Picture 2**, retract the Needle Handle **(2)** fully until it stops and is in position as illustrated in

Picture 2. Unlock the Needle Adjuster (3) by turning counter-clockwise and slide the Adjuster directly below the Needle Handle (2) as shown in **Picture 2**. Turn the Needle Adjuster (3) clockwise to lock the Needle Handle (2)

- ! If any irregularities are noted, the instrument must be disposed of and replaced prior to use.
- ! **Warning: these adjustments must be followed prior to each insertion into the endoscope's working channel, and for each subsequent puncture.**

Handle illustration showing markings, unlocked and locked positions

Markings on the Sliding Sheath Handle (7) and on the Needle Handle (2) indicate if the Adjusters are unlocked to facilitate movement and adjustments, or if they are locked in a safe and locked position.


	
<u>Unlocked</u> Adjuster is positioned on the circle symbol. The unlocking function is performed by turning the Adjuster counter-clockwise.	<u>Locked</u> Adjuster is positioned on the darkened spot. The locking function is performed by turning the Adjuster clockwise.
Picture 3	

Introducing the SonoTip® TopGain EBUS needle into the Endoscope


- When using a Fujifilm bronchoscope, please note before introducing the SonoTip® EBUS needle that the standard enclosed Fujifilm adapter has to be attached as follows:

- Make sure that the Click-Lock Connector **(8)** is opened. (See **Picture 4**).


- Slowly advance the sheath **(9)** with the metal tube which is distally situated on the handle into the endoscope's working channel and Click-Lock **(8)** the handle body onto the scopes working channel. (See **Picture 5**)

- ! **Note:** Ensure that the Needle **(10)** and Nitinol Stylet are securely retracted into the Plastic Sheath **(9)** before advancing the instrument into the accessory channel of the endoscope.
- Adjust the Sheath **(9)** length to the endoscope's working length by unlocking the blue Sheath Adjustor **(6)** and advancing the sliding Sheath Handle **(7)** until the tip of the Sheath **(9)** reaches the desired position. Please ensure that the distal section of the Endoscope is straightened for this purpose.
- ! **Note:** For safe needle advancement, insure that the distal sheath is macroscopically visible on the monitor.
- Secure the sheath length by locking the Sheath Adjustor **(6)**, as shown in **Picture 3**.

Advancing the Needle

- Unlock the Needle Adjustor **(3)** and slowly advance the Needle Handle **(2)** to extend the needle out of the Sheath **(9)** and towards the desired target lesion.
- **Loosen the round stylet from the proximal Luer-Lock of the Needle Handle (2) and retract the stylet not more than 5 mm.**
- Continue to advance the needle **(10)** slowly into the target lesion while maintaining constant visualization via the ultrasonic field. The needle projection length can be set at a maximum projection depth by locking the Needle Adjustor **(3)**. **Picture 6** shows the correct adjustment for a needle projection depth of 2 cm.

- ! Avoid abrupt movements while puncturing and projecting the needle in order not to harm proximate tissue.
- ! If resistance is felt at any time during use, do not exert excessive force as this may result in injury to the patient, endoscope and/or instrument.
- With the needle tip fully embedded into the target lesion, remove the Nitinol Stylet **(1)** completely from the Needle **(10)** and set it aside for reintroduction.
- ! Remove the Stylet carefully to avoid the splashing of possible contaminated tissue/body fluids to minimize risk of infection risk
- Attach the prepared aspiration syringe to the proximal end of the SonoTip® TopGain EBUS and luer lock it securely into place. Release the stopcock on the syringe $\frac{1}{4}$ turn (aligning it parallel with the syringe). This will produce the negative suction pressure needed to aspirate the cell tissue into the needle's tip and lumen.
- Gently advance/retract the Needle **(10)** into the biopsy site to insure adequate sampling using the Needle Handle **(2)**.
- ! **Caution:** Avoid removing the needle **(10)** from the biopsy site during biopsy. This may compromise the integrity of the sample and diagnosis of the biopsy.
- With the needle **(10)** embedded in the biopsy site, close the stopcock on the syringe $\frac{1}{4}$ turn.
- Fully retract the Needle Handle **(2)** so that the Needle **(10)** is removed from the target lesion. Securely house the needle within the Sheath (9) by locking the Needle Adjuster **(3)** in the position as shown in **Picture 2**.
- Open the Click-Lock Connector **(8)** from the endoscope and gently remove the SonoTip® TopGain EBUS from the endoscope.
- Direct the needle tip onto a slide and unlock the syringe plunger to eject the biopsy from the needle **(10)** onto the slide.
- ! The presence of a cytologist during extraluminal endoscopic ultrasound-guided FNB may reduce the number of fine needle biopsies passes required from each lesion.
- If necessary, reinsert the nitinol stylet **(1)** into the needle **(10)** to help eject any remaining aspirate onto the slide. A new syringe with heparinized saline may also be used to flush the needle **(10)** of any excess tissue.
- Once the desired amount of biopsy has been collected, dispose of the entire needle system and syringe according to legal requirements and your facility's infection control protocols.

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Medi-Globe GmbH, Medi-Globe-Str. 1-5, 83101 Achenmühle - Germany
Medi-Globe GmbH