




Save Time. Care Better.

Streamline your practice with the
hassle-free CapsoCam Plus® System.

The CapsoCam Plus System is a sophisticated, user-friendly capsule endoscopy solution that provides a **full 360° image** of the small bowel mucosa - the most comprehensive capsule endoscopy visualization available on the market!

The thorough imaging of the CapsoCam Plus System is enhanced by its wire-free design, which allows for flexibility in patient scheduling, while also providing unprecedented patient freedom and comfort during the examination.

Now, it is the only capsule endoscopy system that is indicated for patients ages 2 and above.



How it Works



CapsoRetrieve[®]

CapsoAccess[®]

CapsoView[®]



1 Patient swallows the vitamin-sized **CapsoCam Plus**.

2 Patient retrieves the **CapsoCam Plus** with the **CapsoRetrieve[®]** tool.

3 Exam data is easily downloaded with the **CapsoAccess[®]** system.

4 Physician quickly plays back and reviews patient data with the **CapsoView[®]** application.



Streamlined Procedure Saves Time and Money

CapsoCam Plus captures a 360° panoramic lateral view of the small bowel using an on-board storage system that eliminates the need for belts, wires, leads and sticky pads.

Removing the need for external receiver equipment, the CapsoCam Plus System comes with zero capital investment. Additionally, administration of a CapsoCam Plus capsule takes less time and doesn't require patients to return bulky equipment that needs charging and maintenance, allowing for more efficient use of staff and exam room time.



Wire-free Design is Patient Preferred¹

CapsoCam Plus patients are completely unburdened by uncomfortable data recorders. During the procedure, patients can go about their daily activities while experiencing unprecedented comfort, convenience, and freedom. In fact, a recent study indicated that **80% of patients prefer CapsoCam¹** over a system with a conventional data recorder.

The patient simply passes the capsule with a normal bowel movement, retrieves it using a specially designed retrieval kit, CapsoRetrieve, and returns it to the clinic by mail or in person.



Data Download is Hassle-Free

The CapsoAccess capsule data access system allows you to seamlessly extract *in vivo* images and other data from the CapsoCam Plus capsule after the procedure. To transfer data from the capsule to your computer, simply place the capsule into the CapsoAccess system, close the lid.



User-friendly Software Makes Reading Faster and Easier

CapsoView software provides a convenient, easy-to-use diagnostic user interface for viewing the CapsoCam Plus images – with several playback options and simple, efficient report generation.



capsovision.com



The CapsoCam Plus video capsule system is intended for visualization of the small bowel mucosa in patients ages 2 and above. It may be used as a tool in the detection of abnormalities of the small bowel.

To learn more or to set up a demonstration, contact your local CapsoCam representative.

Information contained in this leaflet is for use outside the United States only.

¹Patient preference survey administered to more than 100 patients during multi-center North American trial in support of FDA 510(k).
©2017 CapsoVision, Inc. All rights reserved. ADV-2554 Rev C 10/17