

| Bougie Dilator Set

☒ Dilation



Usage

- It is used for treating with esophageal stenosis and cardia stenosis.

Features

- Made by medical macromolecule, with proper rigidity.
- Equipped with radiopaque mark at the tip, easy to observe under X-ray
- Equipped with nitinol super-elastic guide wire, reusable cleaning brush, flusher syringe.

Specifications

Model	Bougie Diameter(mm)	Length(mm)	Characteristics
VDK-BO	5, 7, 9, 11, 13, 15	850	Each set is equipped with 6 bougies, 1 guide wire, 1 reusable cleaning brush and 1 flusher syringe.