


| Disposable Non-electric Jumbo Biopsy Forceps

☒ Biopsy ☒ Polyp Retrieval



Φ3.0MM

Usage

- Tissue sampling in digestive tracts.

Features

- Obtain more adequate biopsy specimens with superior quality.
- Yield better diagnostic accuracy with less tissue acquisition.
- Remove small, sessile colorectal polyps efficiently.

Specifications

Model	O.D.(mm)	Length(mm)	Working Channel(mm)	Characteristics
VDK-FB-30-180-O-O-B1	3.0	1800	≥3.2	without spike, uncoated
VDK-FB-30-180-O-P-B1	3.0	1800	≥3.2	without spike, coated
VDK-FB-30-230-O-O-B1	3.0	2300	≥3.2	without spike, uncoated
VDK-FB-30-230-O-P-B1	3.0	2300	≥3.2	without spike, coated
VDK-FB-30-180-O-P1-B1	3.0	1800	≥3.2	without spike, coated with mark
VDK-FB-30-230-O-P1-B1	3.0	2300	≥3.2	without spike, coated with mark
VDK-FB-30-180-W-P1-B1	3.0	1800	≥3.2	with spike, coated with mark
VDK-FB-30-230-W-P1-B1	3.0	2300	≥3.2	with spike, coated with mark
VDK-FB-30-180-W-O-B1	3.0	1800	≥3.2	with spike, uncoated
VDK-FB-30-230-W-O-B1	3.0	2300	≥3.2	with spike, uncoated
VDK-FB-30-180-W-P-B1	3.0	1800	≥3.2	with spike, uncoated
VDK-FB-30-230-W-P-B1	3.0	2300	≥3.2	with spike, uncoated