


| Gastrointestinal and Biliary Balloon Catheter

☒ Dilation ☒ ERCP



3-STAGE

Usage

- Used for dilation treatment or auxiliary dilation treatment of digestive tract.

Features

- Radiopaque markers---2 X-ray markers on opposite ends to indicate the exact position of the balloon.
- Optimized Passability—New balloon folding technology for a lower profile leading to consistent ease of passage through the scope.
- Varieties of balloon sizes meet all the requirements from different GI &ERCP procedures.

Specifications

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter O.D. (mm)	Catheter Length (mm)	Endoscopic Channel (mm)	Working Pressure (atm)
VDK-TBD-4/5/6-40-90	4/5/6	40	1.8	900	≥2.0	3/6/10
VDK-TBD-4/5/6-40-120	4/5/6	40	1.8	1200	≥2.0	3/6/10
VDK-TBD-4/5/6-40-180	4/5/6	40	1.8	1800	≥2.0	3/6/10
VDK-TBD-4/5/6-40-200	4/5/6	40	1.8	2000	≥2.0	3/6/10
VDK-TBD-6/7/8-40-90	6/7/8	40	2.3	900	≥2.8	3/7/10
VDK-TBD-6/7/8-40-120	6/7/8	40	2.3	1200	≥2.8	3/7/10
VDK-TBD-6/7/8-40-180	6/7/8	40	2.3	1800	≥2.8	3/7/10
VDK-TBD-6/7/8-40-200	6/7/8	40	2.3	2000	≥2.8	3/7/10
VDK-TBD-6/7/8-60-90	6/7/8	60	2.3	900	≥2.8	3/7/10
VDK-TBD-6/7/8-60-120	6/7/8	60	2.3	1200	≥2.8	3/7/10
VDK-TBD-6/7/8-60-200	6/7/8	60	2.3	2000	≥2.8	3/7/10
VDK-TBD-8/9/10-40-90	8/9/10	40	2.3	900	≥2.8	3/5.5/9
VDK-TBD-8/9/10-40-120	8/9/10	40	2.3	1200	≥2.8	3/5.5/9
VDK-TBD-8/9/10-40-200	8/9/10	40	2.3	2000	≥2.8	3/5.5/9
VDK-TBD-8/9/10-60-90	8/9/10	60	2.3	900	≥2.8	3/5.5/9
VDK-TBD-8/9/10-60-120	8/9/10	60	2.3	1200	≥2.8	3/5.5/9
VDK-TBD-8/9/10-60-200	8/9/10	60	2.3	2000	≥2.8	3/5.5/9
VDK-TBD-10/11/12-40-90	10/11/12	40	2.3	900	≥2.8	3/5/8
VDK-TBD-10/11/12-40-120	10/11/12	40	2.3	1200	≥2.8	3/5/8
VDK-TBD-10/11/12-40-200	10/11/12	40	2.3	2000	≥2.8	3/5/8
VDK-TBD-10/11/12-60-90	10/11/12	60	2.3	900	≥2.8	3/5/8
VDK-TBD-10/11/12-60-120	10/11/12	60	2.3	1200	≥2.8	3/5/8
VDK-TBD-10/11/12-60-200	10/11/12	60	2.3	2000	≥2.8	3/5/8
VDK-TBD-12/13.5/15-40-90	12/13.5/15	40	2.3	900	≥2.8	3/5/7
VDK-TBD-12/13.5/15-40-180	12/13.5/15	40	2.3	1800	≥2.8	3/5/7
VDK-TBD-12/13.5/15-40-200	12/13.5/15	40	2.3	2000	≥2.8	3/5/7
VDK-TBD-18/19/20-60-180	18/19/20	60	2.3	1800	≥3.2	3/4.5/6
VDK-TBD-18/19/20-60-200	18/19/20	60	2.3	2000	≥3.2	3/4.5/6
VDK-TBD-18/19/20-60-230	18/19/20	60	2.3	2300	≥3.2	3/4.5/6
VDK-TBD-18/19/20-80-90	18/19/20	80	2.3	900	≥3.2	3/4.5/6
VDK-TBD-18/19/20-80-180	18/19/20	80	2.3	1800	≥3.2	3/4.5/6
VDK-TBD-18/19/20-80-200	18/19/20	80	2.3	2000	≥3.2	3/4.5/6
VDK-TBD-18/19/20-80-230	18/19/20	80	2.3	2300	≥3.2	3/4.5/6