


| Gastrointestinal and Biliary Balloon Catheter

☒ Dilation ☒ ERCP



Usage

- Used for dilation treatment or auxiliary dilation treatment of digestive tract.

Features

- High pressure-resistant materials ensure safe dilation.
- The cylindrical balloon is not easy to deform, and after dilation, it is convenient to observe.
- Various balloon size for choice meets different clinical demand.

Specifications

Dilation balloon for ERCP

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-06-40-200	6	40	2000	5	Suitable for Endoscopic channel $\geq 2.8\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-08-40-200	8	40	2000	5	
VDK-BD-10-40-200	10	40	2000	5	
VDK-BD-12-40-200	12	40	2000	5	
VDK-BD-14-40-200	14	40	2000	5	
VDK-BD-16-40-200	16	40	2000	5	

Dilation balloon for PTCS (percutaneous transhepatic cholechofibroscopy)

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-06-40-90	6	40	900	5	Suitable for Endoscopic channel $\geq 2.0\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-08-40-90	8	40	900	5	
VDK-BD-10-40-90	10	40	900	5	Endoscopic channel $\geq 2.8\text{mm}$

Dilation balloon for Cardia

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-25-80-90	25	80	900	2	Not suitable for Endoscope use
VDK-BD-30-80-90	30	80	900	2	
VDK-BD-35-80-90	35	80	900	2	Guide wire with diameter 0.035" can pass through.
VDK-BD-40-80-90	40	80	900	2	

Dilation balloon for Esophagus, Pylorus and Intestinal tract

Model	Balloon O.D. (mm)	Balloon Length (mm)	Catheter Length (mm)	Working Pressure (atm)	Remark
VDK-BD-10-40-180	10	40	1800	5	Suitable for Endoscopic channel $\geq 2.0\text{mm}$; Guide wire diameter: 0.035"
VDK-BD-12-60-180	12	60	1800	5	
VDK-BD-14-60-180	14	60	1800	5	
VDK-BD-16-80-180	16	80	1800	5	
VDK-BD-18-60-180	18	60	1800	4	
VDK-BD-18-80-180	18	80	1800	4	Endoscopic channel $\geq$ 3.2mm; Not suitable for Guide wire.
VDK-BD-20-60-180	20	60	1800	4	