


| Disposable Hot Biopsy Forceps

☒ ESD

☒ Hemostasis

☒ Polyp Retrieval



Usage

- The product should be equipped with high-frequency instrument, mainly intend to conduct electricity burning, coagulation and hemostasis.

Features

- Easy and comfortable to use.
- With smooth tube, easy to pass through the endoscope channel.
- With puncture-proof insulation, high safety.

Specifications

Model	O.D. (mm)	Length (mm)	Endoscope Channel (mm)	Characteristics
VDK-FD-23-160-A	2.3	1600	≥2.8	Without spike
VDK-FD-23-180-A	2.3	1800	≥2.8	
VDK-FD-23-230-A	2.3	2300	≥2.8	
VDK-FD-23-260-A	2.3	2600	≥2.8	
VDK-FD-23-160-A1	2.3	1600	≥2.8	With spike
VDK-FD-23-180-A1	2.3	1800	≥2.8	
VDK-FD-23-230-A1	2.3	2300	≥2.8	
VDK-FD-23-260-A1	2.3	2600	≥2.8	