


| Stone Extraction Basket

☒ ERCP

☒ Foreign Body Retrieval



***EASY CAPTURE , EASY RELEASE ,
EASY PUSH , OPTIONAL LITHOTRIPSY***

Usage

- It is used for removal of gallstone in biliary duct, or foreign matters in upper and lower digestive tract.

Features

- Innovative handle design, with functions of push, pull and rotation, easier to grasp gallstone and foreign body.
- Convenient to inject contrast medium with injection port on handle.
- Made by advanced alloyed materials, ensures good shape retention even after difficult stone removal.
- Accept customization, can meet different requirements.



Diamond shape



Spiral shape



Semi-circle shape



Tipless shape

Specifications

	Model	O.D. (mm)	Length (mm)	Basket Width (mm)	Endoscope Channel (mm)	Characteristics
	VDK-BAS-18-70-15-S4-A	1.8	700	15	≥ 2.0	Diamond shape
	VDK-BAS-18-70-20-S4-A	1.8	700	20	≥ 2.0	
	VDK-BAS-23-200-25-S4-A	2.3	2000	25	≥ 2.8	
★	VDK-BAS-23-200-30-S4-A	2.3	2000	30	≥ 2.8	
	VDK-BAS-23-230-30-S4-A	2.3	2300	30	≥ 2.8	
	VDK-BAS-23-200-35-S4-A	2.3	2000	35	≥ 2.8	
●	VDK-BAS-18-70-15-N4-B	1.8	700	15	≥ 2.0	Spiral shape
●	VDK-BAS-18-70-20-N4-B	1.8	700	20	≥ 2.0	
	VDK-BAS-23-200-20-N4-B	2.3	2000	20	≥ 2.8	
	VDK-BAS-23-200-25-N4-B	2.3	2000	25	≥ 2.8	
	VDK-BAS-23-200-30-N4-B	2.3	2000	30	≥ 2.8	
	VDK-BAS-18-70-15-N4-C	1.8	700	15	≥ 2.0	Semi-circle shape
	VDK-BAS-18-70-20-N4-C	1.8	700	20	≥ 2.0	
	VDK-BAS-23-200-20-N4-C	2.3	2000	20	≥ 2.8	
★	VDK-BAS-23-200-25-N4-C	2.3	2000	25	≥ 2.8	
	VDK-BAS-23-200-30-N4-C	2.3	2000	30	≥ 2.8	
●	VDK-BAS-18-70-15-N4-D	1.8	700	15	≥ 2.0	Tipless shape
●	VDK-BAS-18-70-20-N4-D	1.8	700	20	≥ 2.0	
	VDK-BAS-23-200-20-N4-D	2.3	2000	20	≥ 2.8	
	VDK-BAS-23-200-25-N4-D	2.3	2000	25	≥ 2.8	

- ★ Common specifications for ERCP
- Common specifications for Cardia